

Pico Green DNA Quantitation

Reagents:

1x TE Procedure calls for “assay buffer.” Use water if diluted samples will be needed for PCR
Prepare from 20X concentrate (stored in -20)
Make ~1mL for each sample you are going to analyze, including duplicates and standards.

Picogreen Prepare 0.5 mL for each sample to be analyzed, including duplicates and standards.
Mix from 200X stock stored in -20. Make dilutions in plastic as Picogreen can sorb to glass. Cover in foil and store in a cool place. Use within 2-3 hours of preparing

2 $\mu\text{g}/\text{mL}$ DNA

Dilute from 50 x stock of λ DNA stored in -20. Amount varies based on standard curve. The table below is for a standard curve set up in a well plate.

DNA ng/mL	2 $\mu\text{g}/\text{mL}$ DNA λ μL	TE μL	Picogreen μL
0	0	100	100
10	1	99	100
50	5	95	100
100	10	90	100
200	20	80	100
500	50	50	100

Sample Prep:

For soil samples prepare a 0.01 dilution of extracted DNA (5 μL in 495 μL water). Add 100 μL of diluted sample to each well and add 100 μL of picogreen dye.

Fluorometer Settings:

Turn on the fluorometer and the computer and allow the fluorometer to warm up for 15 minutes before using.

Under “Instrument” and “Read” the settings should be:

Fluorescence

Excitation 480 Slit Width 15

Emission 520 Slit Width 20

Integrate time 20 seconds

The local contact for the fluorometer is John Shannon jds1c@virginia.edu